DIAGNOSTICS OF PLANT VIRUSES TRAINING MANUAL



Training Workshop on Diagnostics of Plant Viruses (Project No. AGF/CRO/11/007/REG) IPB, UPLB | 17-28 August 2015

aseanet





















DIAGNOSTICS OF PLANT VIRUSES

Marita S. Pinili Lolita M. Dolores Maricel C. Gonzales Raol P. Pamiloza Yron M. Retuta

REMINDER: This manual was prepared/written exclusively for the Training Workshop on Diagnostics of Plant Viruses and may not be reproduced, in whole or in part, in any form or by any means electronic or mechanical, including photocopying, recording, or by any information storage and retrieval system, without the written permission of the ASEAN Plant Health Cooperation Network of the ASEANET (APHCN-ASEANET).

TABLE OF CONTENTS

TABLE OF CONTENTS	Page No.
Sample Collection	1
Indirect ELISA	3
DOT-BLOT ELISA	4
Compound ELISA (BBTV Detection)	5
Total Plant DNA Extraction (Dellaporta Miniprep)	8
CTAB DNA Extraction (For Banana Samples)	9
RNA Extraction	10
Gel Electrophoresis	11
Gel Documentation and Analysis of the Results	12
Polymerase Chain Reaction (PCR)	14
PCR Troubleshooting	15
PCR Cocktail & Thermal Profile Using Degenerate Primers for Begomoviruses	16
PCR Cocktail & Thermal Profile for Banana bunchy top virus (BBTV)	17
PCR Cocktail & Thermal Profile Using Degenerate Primers for Potyviruses	18
PCR Cocktail and Thermal Profile for Papaya ringspot virus (PRSV)	19
Mechanical Inoculation	20
Insect Transmission of Virus	22
Buffer and Solutions	25
Indirect ELISA	25
PBS Buffer (pH7.4)	25
Coating Buffer (pH9.6)	25
Washing Buffer (PBS-T)	25
Blocking Buffer	25
Antibody Buffer	25
Substrate Buffer (pH 9.8)	25
Dot Blot ELISA	26
TBS Buffer (pH8.0)	26
Washing Buffer (TBS-Tween 20)	26
Blocking Buffer	26
Substrate Buffer	26
Nitro Blue Tetrazolium (NBT) Solution	26
5-bromo-4-chloro-3-indolyl phosphate (BCIP)	26
Compound ELISA	27
General Extract Buffer (GEB 1X)	27
Carbonate Coating Buffer (1X)	27
PBST Buffer (Wash Buffer) (1X)	27
ECI Buffer (1X)	27
PNP Buffer (1X)	27

Page

	No
Total DNA Extraction	28
1 M Tris-Base	28
0.5 M EDTA	28
5 M NaCl	28
Dellaporta Extraction Buffer	28
5M Potassium Acetate (KAC)	28
20 % SDS	28
CTAB Extraction Buffer	28
TE Buffer	29
Gel Electrophoresis	29
10x TBE Buffer	29
0.5x TBE Buffer	29
10X Loading Dye	29
6X Loading Dye	29
1Kb plus DNA Ladder (0.1 μg/ml)	29
1.2% Gel	29
GelRed Solution	30
PCR	30
10mM dNTPs	30
Reconstitution of primers	30
Mechanical Inoculation	31
Phosphate Buffer	31
0.01M Phosphate Buffer, pH 7.2	31
Phospahte Buffer (10X)	31
0.01M Phosphate Buffer, pH 7.2	31
Forms	32
ELISA FORM	32
PCR Form for Begomoviruses	33
PCR Form for BBTV	34
RT-PCR Form for Potyviruses	35
RT-PCR Form for PRSV	36
Power Point Presentation/s	37
Detection of Plant Viruses Through Serological Assay (ELISA)	37

SAMPLE COLLECTION

Materials:

- Disposable gloves
- Scissors
- Collection plastics
- Newspaper/Paper towels
- Marker
- Masking tape/Rubber bands
- Cotton
- Tissue paper
- 70% alcohol

Procedure:

1. Cut sample tissue from the plant with clean scissors.

Note: Wipe clean scissors with 70% alcohol in every collection of sample.

2. Put in collection plastic with wet cotton. Tie a knot to close or use a rubber band. Another way is to wrap the sample in newspaper or paper towel then secure using a masking tape.

Note: Properly label each sample. Don't forget to include the place and date of collection.

- 3. Fill-up sampling details in the collection form.
- 4. Extract the DNA or RNA immediately after collection or store samples in a refrigerator then process it the next day.
 - Note: RNA can be easily degraded, that is why it is important to immediately process the samples. Also some plant tissues oxidized quickly like banana leaves. However, some plant samples can be stored in silica gel like tomato, squash and pepper leaves.

		Σ												
	Date of Collection	Feb 10, 201												
	Estimated %Infection	30												
	Symptoms	Leaf curling, mosaic												
PLE\$	Farm Size	500 <u>sq.m</u> .												
ANT \$AM	Age of Plant (days)	120												
LIST OF PL	Variety	Suprema												
	Farmer's or Farm's Name	DA Paraiso												
	Place of Collection	Tayug, Pangasinan												
	Crop	Squash												
	Sample No.	-												

INDIRECT ELISA

Materials:

- Blade or scissors
- Collection plastic (4x6in)
- Paper towels
- Wash Bottle
- Pipettes and Tips
- Disposable gloves

Equipment:

- ELISA Reader
- Weighing balance
- Incubator (Set at 37 °C)

Buffers and Solutions:

- PBS Buffer (pH7.4)
- Coating Buffer (pH9.6)
- Blocking Buffer
- Antibody Buffer
- Washing Buffer (PBS-T)
- Substrate Buffer (pH 9.8)
- Virus-specific antibody
- Goat anti-rabbit enzyme conjugate (GARAP)
- p-ntrophenylphosphatae (pNPP)
- 3M NaOH (Stop solution)

- 1. Grind leaf samples in coating buffer in a dilution of 1:10.
- 2. Load 200 µl plant sap per well of the ELISA plate.
- 3. Incubate overnight at 4 °C.
- 4. After incubation, remove the plant sap from the wells. Fill the wells with washing buffer; let it stand for 5 min then empty the plates. Repeat 3 times.
- 5. Add 300 μI blocking solution to each well and incubate at room temperature for 1 hr.
- 6. Repeat step 4 (Washing).
- 7. Add 100 μ l virus-specific antibody with a dilution of 1:200 in antibody buffer to each well.
- 8. Incubate at 37 °C for 2-3 hrs.
- 9. Repeat washing as in step 4.
- 10. Add 100 µl goat anti-rabbit enzyme conjugate (GARAP) with a dilution of 1:1000 in antibody buffer.
- 11. Incubate at 37 °C for 2-3 hrs.
- 12. Repeat washing step.
- 13. Add 100 μl of p-nitrophenylphosphatase (pNPP) in substrate buffer to each well.
- 14. Incubate at room temperature for 30-60 min. Observe color reaction.
- 15. Stop reaction by adding 50 µl 3M NaOH.
- 16. Assess results by:
 - Visual observation
 - Absorbance Reading at 405 nm using an ELISA Reader.

DOT-BLOT ELISA

Materials:

- Blade or scissors
- Collection plastic (4x6in)
- Paper towels
- Wash Bottle
- Pipettes and Tips
- Disposable gloves

Equipment:

• Weighing balance

Buffers and Solutions:

- Nitrocellulose membrane (45 µm)
- TBS Buffer (pH8.0)
- Blocking Buffer
- Washing Buffer (TBS-Tween 20)
- Substrate Buffer (pH 8.0)
- Virus-specific antibody
- Goat anti-rabbit enzyme conjugate (GARAP)
- Nitro Blue Tetrazolium (NBT)
- 5-bromo-4cholro-3-indolyl phosphate (BCIP)
- 1.5% Sodium hypochlorite (Stop solution)

- 1. Gently draw a grid pattern on the nitrocellulose membrane (NCM) with a pencil and ruler. Dip membrane onto PBS buffer then air-dry before using.
- 2. Grind samples at a ratio of 1:20 (w/v) in TBS buffer.
- 3. Drop 2-3 µl plant sap onto the membrane. Incubate overnight at room temperature.
- 4. Evenly cover the membrane with blocking buffer in a plastic container with cover then incubate at room temperature for 1 hr.
- 5. Remove blocking buffer. Dry container with clean paper towel.
- 6. Add 1:200 (v/v) of virus specific antibody in TBS buffer.
- 7. Incubate at room temperature for 1-2 hr.
- 8. After incubation, remove antibody and buffer then rinse 3 times with TBS-T with 10 min interval.
- Add 1:1000 (v/v) goat anti-rabbit enzyme conjugate (GARAP) in blocking buffer.
- 10. Incubate at room temperature for 1-2 hr.
- 11. Remove conjugate buffer. Incubate in TBS buffer at room temperature for 30 min.
- 12. Remove buffer. Soak membrane in 5 ml substrate buffer with 33 μl NBT and 16.5 μl BCIP.
- 13. Observe color reaction. After color reaction (purple or blue) developed in the positive check, pour off the substrate buffer.
- 14. Add 1.5 % sodium hypochlorite and incubate at room temperature for 5 min to eliminate color background and stop the reaction.

COMPOUND ELISA

(BBTV Detection)

Materials:

- Blade or scissors
- Collection plastic (4x6in)
- Paper towels
- Wash Bottle
- Pipettes and Tips
- Disposable gloves

Equipment:

- ELISA Reader
- Weighing balance

Procedure:

1. Prepare capture antibody.

Note: All antibodies and enzyme conjugates should be prepared in a container made of a material such as polyethylene or glass that does not readily bind antibodies. Do not use polysterene, polypropylene or polycarbonate

Prepare the volume of carbonate coating buffer needed for the test. *Example: If the dilution given on the bottle of concentrated capture antibody is 1:200, and you are preparing 10 ml of the capture antibody solution, you should mix 10 ml of carbonate coating buffer with 50 µl of concentrated capture antibody. Mix the prepared antibody solution thoroughly and use immediately.*

- 2. Coat plate. Pipette 100 μl of the prepared capture antibody solution into each well.
- 3. Incubate plate. Cover the plate with cling plastic wrap and incubate in a humid box for 4 hrs at room temperature or overnight in the refrigerator (4°C). *Do not store coated plates longer than 24 hrs.*
- 4. Wash plate. Empty the wells into a sink or container. Fill the test wells completely with 1X PBST, and then quickly empty them again. Repeat 2 more times.

Hold the plate upside down and tap firmly on folded paper towel to remove excess liquid or if you have an automated plate washer, calibrate the washer and wash the plate 3x.

Buffers and Solutions:

- General Extract Buffer (GEB 1X)
- Carbonate Coating Buffer (1X)
- PBST Buffer (Wash Buffer) (1X)
- ECI Buffer (1X)
- PNP Buffer (1X)

5. Grind and dilute the samples. Select samples showing symptoms. Young leaf tissue is recommended. Seed, stem and other tissue can also be tested depending on the crop.

Grind plant tissue in General Extract Buffer (GEB) at 1:10 ratio (*tissue weight in g: buffer volume in ml*). You can use mortar and pestle or any grinding devices. Be sure to wash and rinse the grinding device thoroughly between samples.

- Dispense the sample. Following the loading diagram (see sample diagram), dispense 100 µl of prepared sample into sample wells. Dispense 100 µl of positive control into positive control wells, and dispense 100 µl of sample extraction buffer into buffer wells.
- 7. Incubate plate. Set the plate, cover it with cling wrap and incubate in humid box for 2 hrs at room temperature or overnight in the refrigerator (4°C).
- 8. Prepare enzyme conjugate.

Note: always prepare enzyme conjugate within 10 min before use. Bottles of alkaline phosphatase enzyme conjugate and detection antibody are supplied as a concentrate and must be diluted with ECI buffer before use. Please follow the recommended dilution stated in the product.

Example: If the dilution given on the bottles of concentrated detection antibody (A) and alkaline phosphatase enzyme (B) conjugate is 1:200, and you are preparing 10 ml of enzyme conjugate solution, you should first dispense 10 ml of ECI buffer. Then, add 50 μ l of A and 50 μ l of B to the ECI buffer.

After adding the reagents from A and B, mix thoroughly the enzyme conjugate solution well.

- 9. Wash plate. When sample incubate is complete, wash the plate. Use the quick flipping motion to dump the wells into a sink without mixing the contents. Fill all the test wells completely with 1X PBST, and then quickly empty them. Repeat for 2 more times. After washing tap the plate firmly on folded paper towel to remove all droplets of wash buffer.
 - Note: It is important to inspect the test wells. It should be free from any plant tissue and dirt. If plant tissue is present repeat the wash step and tap firmly. Avoid touching the bottom of the plate.
- 10. Add enzyme conjugate. Dispense 100 µl of prepared enzyme conjugate solution per well.
- 11. Incubate the plate. Cover the plate with cling wrap and incubate in humid box for 2 hrs at room temperature.

12. Prepare PNP solution. Each PNP tablet will make 5 ml of PNP solution, at a concentration of 1 mg/ml about enough for five 8-well strips. About 15 min before the end of incubation step, measure 5 ml of room temperature 1X PNP buffer for each tablet you will be using. Then without touching the tablets, add the PNP tablets to the buffer.

Note: Do not touch the PNP tablets or expose the PNP solution to strong light. Light or contamination could cause background color in negative wells.

- 13. Wash plate. Follow the same wash procedure previously described.
- 14. Add PNP substrate. Dispense 100 µl of PNP substrate into each testwell.
- Incubate plate. Cover the plate with cling wrap and incubate in humid box for 60 min at room temperature. Plates should be protected from direct or intense light.
- 16. Evaluate results.

Qualitative interpretation. Examine each well by eye, or measure on a plate reader at 405 nm. Air bubbles should be removed; it can alter results at the time of reading.

Wells in which color develops indicate positive results. Wells in which there is no significant color development indicate negative result. Test results are valid only if positive control wells give a positive result and buffer wells remain colorless.

Results may be interpreted after more than 60 min of incubation as long as negative wells remain virtually clear.

Quantitative interpretation. Appropriate controls should be included for reference, since ELISA values may differ in different microtiter plates due to possible plate-to-plate variation in sensitivity. Overlapping range of specific and non-specific reaction values causes difficulty in interpretation. In this case, it is necessary to include large number of known healthy control samples and determine statistically a threshold level for infection. *To establish thresholds, several authors have used the mean value for healthy controls plus three times their standard deviation (\bar{x} + 3SD).*

Alternatively, values more than twice those healthy controls have been considered infected.

Training Workshop on Diagnostics of Plant Viruses

IPB, University of the Philippines Los Baños | 17-28 August 2015

TOTAL PLANT DNA EXTRACTION

(Dellaporta Miniprep)

Materials:

- Blade or scissors
- Mortar and pestle (sterilized)
- 1.5 mL Eppendorf tubes
- Pipettes and Tips
- Face mask
- Disposable gloves
- Paper towels

Equipment:

Buffers and Solutions:

- Nitrocellulose membrane (45 µm)
- Dellaporta extraction buffer
- 20% Sodium Dodecyl Sulfate (SDS)
- 5M Potassium acetate (KAC)
- Isopropanol
- 80% ethanol
- Sterile distilled water or TE buffer

- Fume hood
- Water bath or Dri-bath incubator
- Centrifuge
- Vortex
- Weighing balance

- 1. Collect 2 leaf discs about 0.9cm in diameter.
- 2. Homogenize in 500 µl Dellaporta extraction buffer using a mortar and pestle.
- 3. Transfer in a 1.5 ml microfuge tube then add 33 μl 20% SDS. Gently invert tube to mix the solution.
- 4. Incubate at 65°C for 10min.
- 5. Add 160 µl 5M KAC and mix by gently inverting the tube.
- 6. Centrifuge at 13,000 rpm for 10 min.
- 7. Transfer supernatant into a new 1.5 ml microfuge tube, avoiding the plant tissue debris.
- 8. Repeat centrifugation. Transfer 500µl supernatant into a new tube.
- 9. Add 0.5 volume (250 µl) isopropanol and invert tube gently.
- 10. Centrifuge at 13,000 rpm for 10min.
- 11. Carefully discard supernatant Note: Make sure that the pellet does not become aspirated. It may be necessary to leave some (about 15 μl) supernatant behind. The following step will delete any problems it may cause.
- 12. Add 500 µl 80% ethanol and centrifuge at 13,000 rpm for 5min.
- 13. Carefully discard as much supernatant as possible.
- 14. Air-dry for 1hr or speed-vacuum for 5min.
- 15. Resuspend pellet in 500 µl sterile distilled water or DEPC-treated water.

CTAB DNA EXTRACTION

(For Banana Samples)

Materials:

- Blade or scissors
- Mortar and pestle (sterilized)
- 1.5 mL Eppendorf tubes
- Pipettes and Tips
- Face mask
- Disposable gloves
- Paper towels

Equipment:

Buffers and Solutions:

- CTAB extraction buffer
- 2-mercaptoethanol
- Phenol:Chloroform:Isoamyl (PCI) 25:24:1
- 95% ethanol
- 70% ethanol
- Sterile distilled water or TE buffer

- Fume hood
- Water bath or Dri-bath incubator
- Centrifuge
- Vortex
- Weighing balance

- 1. Clean working station and pipettors with 70% alcohol.
- 2. Cut 0.5 g leaf sample using sterile blade or clean scissors.
- 3. Homogenize leaf sample using a sterile mortar and pestle with 4 mL CTAB extraction buffer and 8 µl 2-mercaptoethanol.
- 4. Transfer 500 µL leaf extract in a sterile 1.5 mL sterile microfuge tube.
- 5. Incubate in a water bath at 60 °C for 1 hr. Gently swirl tubes every 15 min.
- After incubation, add 400 μL Phenol:Chloroform:isoamyl (PCI 25:24:1). Mix by vortexing for 15 sec.
- 7. Centrifuge at 10,000 rpm for 5 min.
- 8. Collect 500 μ L of the upper phase using a wide-bore tip and transfer it to a new 1.5 mL microfuge tube.
- 9. Precipitate the DNA by adding 2 volumes of ice-cold 95% ethanol. Mix well by gently inverting the tube. Incubate in the freezer for 1 hr or overnight.
- 10. Centrifuge at 10,000 rpm for 15 min.
- 11. Discard the supernatant. Wash the pellet with 500 μ L ice-cold 70% ethanol and gently invert the tube for 3-5 min.
- 12. Short spin then discard the ethanol. Air-dry the pellet for 2-3 min at room temperature in the laminar flow or until the wall of the tube is already dry.
- 13. Resuspend the pellet in 100 uL TE buffer or DNAse-free water.
- 14. Store at -20°C.

RNA EXTRACTION

Materials:

- Blade or scissors
- Mortar and pestle (sterilized)
- 1.5 mL Eppendorf tubes
- Pipettes and Tips
- Face mask
- Disposable gloves
- Paper towels

Equipment:

Buffers and Solutions:

- Trizol® Reagent (Invitrogen)
- Chloroform
- Isopropanol
- 75% ethanol
- RNAse-free water or TE buffer
- RNAse Away
- 70% alcohol

- Fume hood
- Water bath or Dri-bath incubator
- Refrigerated centrifuge
- Vortex

- 1. Clean working station and pipettors with 70% alcohol or RNAse Away.
- 2. Homogenize about 100 mg leaf sample in 1 ml Trizol® Reagent using an icecold mortar and pestle.
- 3. Incubate at room temperature for 5 min.
- 4. Centrifuge at 12,000xg for 10 min at 8°C.
- 5. Transfer the supernatant using a wide-bore tip into a new 1.5 ml microfuge tube then add 200 μl chloroform.
- 6. Shake vigorously by hand for 15 sec.
- 7. Incubate at room temperature for 3 min then centrifuge at 12,000xg for 10 min at 8°C.
- 8. Collect the aqueous phase into a new 1.5 ml microfuge tube then add 500 μl ice-cold isopropanol.
- 9. Incubate at -20 °C for 2 hrs or overnight if precipitate fails to form.
- 10. Discard supernatant and wash pellet with ice-cold 75 % ethanol.
- 11. Centrifuge at 7,500xg for 5 min at 8°C.
- 12. Discard ethanol then air-dry pellet for 10 min in a laminar flowhood.
- 13. Resuspend pellet in 50 µl RNAse-free water or TE buffer.
- 14. Incubate at 55 °C for 10 min before storing at -20 °C or in a biofreezer (-80 °C).

GEL ELECTROPHORESIS

Materials:

- Casting tray
- Gel comb
- Pipettes and Tips
- Disposable gloves

Equipment:

- Gel electrophoresis system
- Weighing balance
- Microwave oven

Procedure:

- 1. Prepare 1.2% agarose gel in 0.5X TBE buffer.
- 2. Completely dissolve the agarose in the buffer using a microwave.
- 3. Let the solution to slightly cool down (~5 min).
- 4. Pour the solution slowly into the casting tray with the gel comb in place. Avoid forming any bubbles.
- 5. Let the agarose gel to solidify (~30 min) then carefully remove the gel comb.
- 6. Place the solidified agarose gel into the electrophoresis unit. Fill in with 0.5X TBE buffer until the gel is fully submerged.
- Pipette several 2 µl loading dye (depending on the number of samples) in a piece of parafilm.
- 8. Mix 5 µI PCR product into the dye by carefully pipetting the solution in and out of the tip. Avoid forming of bubbles.
- 9. Load the sample mixture into each well of the gel. Avoid spilling over the sides of the wells to prevent contamination. *Note: Start loading samples on the second well.*
- 10. Lastly, load 2 µl DNA ladder on the first well.
- 11. Run samples for 30-45 min or until the dye has migrated about 75-80% of the gel.
- 12. After the run, stain the gel by submerging it in the GelRed solution for about 15-30 min.
- 13. Visualize the amplified DNA bands using a gel documentation system.

Buffers and Solutions

- 0.5X TBE
- Agarose
- GelRed
- 1 KB plus DNA Ladder

GEL DOCUMENTATON AND ANALYSIS OF THE RESULTS

Materials:

Equipment:

- Disposable gloves
 - Kimwipes
- Paper towel or tissue

- GelDoc Computer

- 1. After staining the gel, the PCR result can be viewed using a gel documentation system or GelDoc.
- 2. Analyze the result using the DNA ladder.

Kb Plus DNA Ladder	1 Kb DNA Ladder							
bp	bp		bp					
-12,000	12,216 -		-11,198					
	8,144 -		-9.162					
-5,000	6,108 -		= 7,120 = 000					
	4,072 -		- 5,090					
-2.000	3,054 -							
-1,650	10000000							
	2,036 -							
1 000	1,636 -	-						
- 850								
- 650								
- 500								
-400								
- 300	506 -							
- 200								
- 100								
	1							
0.9 µg/lane	0.	9 µg/la	ne					

Example of DNA ladders (Invitrogen)

- M 1 2 3 4 5 6 7 12,000bp 1650bp 1000bp 100bp
- 3. If the PCR run is successful, the PCR product is within the expected DNA sequence size.

Agarose gel showing the expected band size of around 1.5kb after PCR amplification.

4. The presence of the expected band, confirms the presence of the virus.

Note: It is important that in every gel electrophoresis for PCR analysis, **DO NOT FORGET TO INCLUDE THE DNA LADDER.**

POLYMERASE CHAIN REACTION (PCR)

Materials:

- 0.2ml PCR tubes
- 1.5 ml microfuge tubes
- 70% ethanol
- Pipettes and Tips
- Disposable gloves
- Paper towel or tissue

Equipment:

- PCR machine
- Centrifuge

Buffers and Solutions

- 10X PCR buffer
- 50 mM MgCl₂
- 10 mM dNTPs
- 10 µM forward and reverse primers
- RNAse-free water
- Taq Polymerase
- Template DNA

- 1. Before starting:
 - a. Calculate the amount of the PCR reagents needed according to the number of samples to be tested. Give allowance for pipetting error.
 - b. Clean working bench, pipettors and microfuge rack with 70% ethanol.
 - c. Thaw reagents (PCR buffer, MgCl₂, dNTPs, forward and reverse primers, RNAse-free water and DNA samples) on ice.
 - d. Flick tubes to mix its content then short spin.
 - e. Prepare and label 0.2ml PCR tubes.
- 2. Prepare a cocktail mix. Add all the PCR reagents in appropriate amounts in a 1.5 ml microfuge tube.
- 3. Lastly add the Taq polymerase to the cocktail.
- 4. Gently flick the tube to mix the reagents then short spin the tube using a centrifuge.
- 5. Dispense the PCR mix on each labeled 0.2 ml PCR tube. Add the DNA template then gently mix the contents.
- 6. Place tubes in the thermal cycler then begin run with the appropriate PCR program.

PCR Troubleshooting

CAUSE	SOLUTION				
Problem: Faint bands or no PCR	product				
a. Too little DNA template in the reaction	Increase the amount of the DNA template.				
b. Damaged or degraded DNA template	Assure the purity and integrity of the DNA template. Be careful in handling DNA. Prevent freeze-thawing of the DNA samples by preparing aliquots.				
c. Insufficient <i>Taq</i> Polymerase	Increase the DNA polymerase concentration in increments of 0.5 units per 100 µL of reaction.				
d. Insufficient number of cycles	Increase cycle number by 5 to 10 cycles.				
e. Presence of PCR inhibitors	Re-purify DNA samples.				
f. To low MgCl ₂ concentration	Increase magnesium concentration in increments of 0.1mM.				
g. Too long or too short denaturation time.	Adjust denaturation time in increments of 5sec.				
h. Too high annealing temperature	Lower annealing temperature in increments of 2° C. Compute the melting temperature (T _m) of the primers. The annealing temperature should be 5° C less to the primer T _m .				
i. The primer extension period is too short.	Increase extension time in increments of 1 minute.				
Multiple bands or smearing					
a. Too much DNA template	Decrease the amount of template DNA				
b. Too low annealing temperature	Increase annealing temperature in increments of 2°C.				
c. Too high concentration of <i>Taq</i> Pol	Decrease enzyme concentration in increments of 0.5 units per 100-µL reaction.				
d. Magnesium concentration is too high.	Decrease the magnesium concentration in increments of 0.1mM.				
e. Denaturation time is too short or too low	Increase the denaturation time in increments of 5sec. and temperature by 1°C.				
f. Too many cycles	Reduce the cycle number by 5 to 10 cycle.				
g. Extension time is too long	Reduce the extension time in increments of 1min.				
h. Review primer design and composition.	Design new primers.				

PCR Cocktail & Thermal Profile Using Degenerate Primers for Begomoviruses

1. PCR Cocktail

PCR Cocktail Preparation							
Components	Stock Concentration	Final Concentration	1x				
DEPC-treated Water			17.65				
PCR Buffer	10x	1x	2.5				
MgCl ₂	50mM	2.5 mM	1.25				
dNTPS	10mM	0.2 mM	0.5				
Primer F	10µM	0.5 µM	0.5				
Primer R	10µM	0.5 µM	0.5				
Таq	5Ú/µl	0.06	0.1				
DNA	·		2				
Total Reaction Volume			25 µl				
Aliquot 23µl to each DNA tube							

2. PCR Profile

Steps	Temperature (°C)	Time	Number of Cycles
Denaturation	94	1 min	
Annealing	57	2 min	30x
Elongation	72	2 min	
Final Extension	72	10 min	1
Hold	16	∞	-

3. Primer Sequence (amplifies the fragment of DNA-A including the 5' end of CI, IR, V2 and the 5' end of the coat protein (CP), Tsai etal., 2011)

Primer Name	Sequence	Amplicon Size (bp)
PAL1v1978B	GCATCTGCAGGCCCACATBGTYTTHCCNGT	
PAR1c715H	GATTTCTGCAGTTDATRTTHTCRTCCATCCA	~1.5 KD

PCR Cocktail & Thermal Profile for Banana bunchy top virus (BBTV)

1. PCR Cocktail

PCR Cocktail Preparation							
Components	Stock Concentration	Final Concentration	1x				
DEPC-treated Water			17.65				
PCR Buffer	10x	1x	2.5				
MgCl ₂	50mM	2.5 mM	1.25				
dNTPs	10mM	0.2 mM	0.5				
Primer F	10µM	0.5 µM	0.5				
Primer R	10µM	0.5 µM	0.5				
Таq	5Ú/µl	0.06	0.1				
DNA	·		2				
Total Reaction Volume			25 µl				
Aliquot 23µl to each DNA tube							

2. PCR Profile

Steps	Temperature (°C)	Time	Number of Cycles
Initial Denaturation	94	4min	1
Denaturation	94	1 min	I
Annealing	61	1 min	29x
Elongation	72	2 min	
Final Extension	72	10 min	1
Hold	16	∞	-

3. Primer Sequence (amplifies the BBTV-R genome)

Primer Name	Sequence	Amplicon Size (bp)
D11	GGAAGAAGCCTCTCATCTGCTTCAGACARC	
D12	TTCCCAGGCGCACACCTTGAGAAACGAAAG	~1.1 KD

RT-PCR/PCR Cocktail & Thermal Profile Using Degenerate Primers for *Potyviruses*

1. RT-PCR Cocktail

PCR Cocktail Preparation							
Components	Stock Concentration	Final Concentration	1X				
DEPC-treated Water			2.2				
Reaction Mix	2X	1X	5				
Primer F	10µM	0.5 µM	0.2				
Primer R	10µM	0.5 µM	0.2				
SuperScript III RT	-	-	0.4				
RNA			2				
Total Reaction Volume			10 µl				
Aliquot 8 µl to each DNA tube							

2. PCR Profile

Steps	Temperature (°C)	Time	Number of Cycles
cDNA Synthesis	55	30 min	1
Initial Denaturation	94	2 min	1
Denaturation	94	15 sec	
Annealing	55	30 sec	40x
Elongation	68	1 min	
Final Extension	68	5 min	1
Hold	16	∞	-

3. Primer Sequence (amplifies the 3' terminal portion of the genomes of various *potyviruses*, Gibbs and Mackenzie, 1997)

Primer Name	Sequence	Amplicon Size (bp)
Potyvirid 1	CACGGATCCCGGG(T)17VGC	1.6 Kb
Potyvirid 2	ACCACAGGATCCGGBAAYAAYAGYGGDCARCC	~1.0 ND

RT-PCR/PCR Cocktail and Thermal Profile for Papaya ringspot virus (PRSV)

1. RT-PCR Cocktail

PCR Cocktail Preparation			
Components	Stock Concentration	Final Concentration	1X
DEPC-treated Water			2.2
Reaction Mix	2X	1X	5
Primer F	10µM	0.5 µM	0.2
Primer R	10µM	0.5 µM	0.2
SuperScript III RT	-	-	0.4
RNA			2
Total Reaction Volume			10 µl
Aliquot 8 µl to each DNA tube			

2. PCR Profile

Steps	Temperature (°C)	Time	Number of Cycles
cDNA Synthesis	55	30 min	1
Initial Denaturation	94	2 min	1
Denaturation	94	15 sec	
Annealing	58	30 sec	40x
Elongation	68	1 min	
Final Extension	68	5 min	1
Hold	16	∞	-

3. Primer Sequence (amplifies the PRSV-CP gene, Bateson etal., 1994)

Primer Name	Sequence	Amplicon Size (bp)
MB11	GGATCCATGTCCAAAAATGAAGCTGTGGATGCT	000 h a
MB12	TCAATTGGCGCATACCCAGGAGAGT	~900 pb

MECHANICAL INOCULATION

Materials:

- Mortar and pestle
- Wash bottle
- Carborundum (500mesh)/Celite
- Disposable gloves
- Labels and water resistant marker
- Detergent soap

Buffers and Solutions

- 0.01M Phosphate Buffer pH 7.2
- Sodium sulphite (2%)

Procedure:

- 1. Arrange plants to be inoculated.
- 2. Homogenize infected leaf sample in 0.01M Phosphate buffer at 1:10 dilution using a mortar and pestle. Add 2% Sodium sulphite.
- 3. Add a small amount of celite (0.5-1% w/v) onto the inoculum or spread small amount of carborundum into the leaf to be inoculated.
- 4. Moisten two fingers with the inoculum then gently rub onto the first and second fully expanded leaves while supporting it with the other hand.
- 5. Rinse inoculated plants with tap water within 2-5 min.
- 6. Observe test plants for symptoms at least twice a week for a month.
- 7. Label with date of inoculation and name of virus inoculation.

Note: To prevent cross contamination, change or wash thoroughly gloves between samples. Also separate test plants of different samples.

Test plant species	Number of leaves	Remarks
Abelmoschus esculentus	1-2	
Capsicum annuum	2-3	
Chenopodium amaranticolor	3-4	for local symptoms only
Chenopodium quinoa	3-4	for local symptoms only
Cucumis sativus	2 cotyledons	remove leaves, except top leaf
Cucurbita maxima	2 cotyledons	
Datura metel	2-3	
Gomphrena globosa	about 6	
Gossypium herbaceum	1-2	
Nicotiana benthamiana	3-4	
Nicotaina glutinosa	3-4	
Nicotiana tabacum 'Xanthii'	1-2	
Solanum lycopersicum 'Money- maker'	1-2	
Solanum melongena	1-2	
Vigna unguiculata	2 cotyledon	

Optimum stage of most common test plant species.

Source: <u>www.q-bank.eu</u>

INSECT TRANSMISSION OF VIRUS

I. Non - Persistent Mode of Transmission

Materials:

- Insect vector, Aphis gossypii
- Test plants
- Symptomatic plants/Inoculum
- Camel hair brush
- Close container with screen window

- 1. Starve the non-viruliferous insect vectors for 30 minutes before virus acquisition.
- 2. Place the starved insect vectors to symptomatic host plant for 15-20 minutes to acquire virus.
- 3. Transfer the viruliferous insect vectors (15-20 aphids) to the uninfected/clean test plant and allow the insects to feed and transmit the virus for 1-2 hours.
- 4. Manually remove or eliminate the insects by spraying insecticide.
- 5. Transfer the test plants to insect-free/insect-proof cages.
- 6. After 1-2 months, check the test plants for virus infection symptom/s and collect samples for ELISA and/or PCR tests.

INSECT TRANSMISSION OF VIRUS

IIa. Persistent Mode of Transmission (Banana bunchy top virus)

Materials:

- Insect vector, Pentalonia nigronervosa
- Test plant, 1-mo old banana seedling
- Symptomatic plants/Inoculum
- Camel hair brush
- Close container with screen window

Procedure:

- 1. Starve the non-viruliferous insect vectors for 30 minutes before virus acquisition.
- 2. Place the starved insect vectors to symptomatic host plant for 24 hours to acquire virus.
- 3. Transfer the viruliferous insect vectors (15-20 aphids) to the uninfected/clean test plant and allow the insects to feed and transmit the virus for 24 hours.
- 4. Manually remove or eliminate the insects by spraying insecticide.
- 5. Transfer the test plants to insect-free/insect-proof cages.
- 6. After 1-2 months, check the test plants for virus infection symptom/s and collect samples for ELISA and/or PCR tests.

Note: In the case of Banana bunchy top virus, symptom expression may take (up to 3 months) depending on the susceptibility of the host plant.

INSECT TRANSMISSION OF VIRUS

IIb. Persistent Circulative Mode of Transmission (Tomato leaf curl virus)

Materials:

- Aviruliferous whiteflies, *Bemisia* tabaci
- ToLCV infected tomato (source of inoculum)
- Tomato seedlings (healthy)
- Aspirator
- Screen cages

- 1. Collect/aspirate whiteflies and allow having access on ToLCV infected tomato plant for 24-48 hours.
- 2. Remove whiteflies from the source and transfer 10-20 whiteflies to healthy tomato seedling inside insect proof screen cages.
- 3. Allow the viruliferous whiteflies to feed on test plants for 48 hours then remove or spray insects with insecticide.
- 4. Observe the inoculated plants for virus symptoms (leaf curl) after 3-4 weeks until 2 months and collect samples for PCR.

BUFFERS AND SOLUTIONS

Indirect ELISA

PBS Buffer (pH7.4) Sodium Chloride (NaCl) Monobasic Potassium Phosphate (KH ₂ PO ₄) di-Sodium hydrogen Orthophosphate (Na ₂ HPO ₄ x12H ₂ O) Potassium Chloride (KCl) Distilled water	8 g 0.2 g 1.44 g 0.2 g 1000 ml
<u>Coating Buffer (pH9.6)</u> Sodium Carbonate (Na ₂ CO ₃) Sodium Hydrogen Carbonate (NaHCO ₃) Distilled water	1.5 g 2.93 g 1000 ml
<u>Washing Buffer (PBS-T)</u> PBS Buffer Tween 20	1000 ml 0.05%
<u>Blocking Buffer</u> (must be freshly prepared) PBS Buffer BSA or skim milk	500 ml 1%
<u>Antibody Buffer</u> (must be freshly prepared) PBS Buffer Egg albumin	100 ml 0.2%
<u>Substrate Buffer (pH 9.8)</u> Diethanolamine Distilled water Adjust pH to 9.8 then volume to 1000 ml.	97 ml 800 ml

Training Workshop on Diagnostics of Plant Viruses

IPB, University of the Philippines Los Baños | 17-28 August 2015

Dot Blot ELISA

TBS Buffer (pH8.0)	
Tris Distilled water Adjust pH to 8.0	6.057 g 800 ml
NaCl Volume to 1000 ml.	8.766 g
Washing Buffer (TBS-Tween 20) TBS Buffer Tween 20	1000 ml 0.5 ml
Blocking Buffer Skim Milk Glycine TBS-Tween 20	3 g 2 g 100 ml
Substrate Buffer (pH 8.0) Tris Dissolve in 80 ml distilled water then adjust pH to 8.0. NaCl MgCl ₂ x6H ₂ O Volume up to 100 ml.	1.214 g 0.5844 g 0.102 g
Nitro Blue Tetrazolium (NBT) Solution Nitro Blue tetrazolium chloride Dimethylformamide (70%) Store at 4°C.	0.5 g 10 ml
<u>5-bromo-4-chloro-3-indolyl phosphate (BCIP)</u> Disodium salt BCIP Dimethylformamide (100%)	0.5 g 10 ml

Compound ELISA

<u>General Extract Buffer (GEB 1X)</u> Dissolve in 1000 ml of 1X PBST: Sodium sulfite (anhydrous) Polyvinylpyrrolidone (PVP) MW 24-40,000 Sodium azide Powdered egg (chicken) albumin, Grade II Tween-20 Adjust pH to 7.4. Store at 4°C	1.3g 20g 0.2g 2g 20g
Carbonate Coating Buffer (1X) Dissolve in distilled water to 1000 ml Sodium carbonate (anhydrous) Sodium bicarbonate Sodium azide Adjust pH to 9.6. Store at 4°C	1.59g 2.93g 0.2g
PBST Buffer (Wash Buffer) (1X) Dissolve in distilled water to 1000 ml: Sodium chloride Sodium phosphate, dibasic (anhydrous) Potassium phosphate, monobasic (anhydrous) Potassium chloride Tween-20 Adjust pH to 7.4	8g 1.15g 0.2g 0.2g 0.5g
ECI Buffer (1X) Add to 1000 ml 1X PBST: Bovine serum albumin (BSA) Polyvinylpyrrolidone (PVP) MW 24-40,000 Sodium azide Adjust pH to 7.4 Store at 4°C	2g 20g 0.2g
<u>PNP Buffer (1X)</u> Dissolve in 800 ml distilled water: Magnesium chloride hexahydrate Sodium azide Diethanolamine Adjust pH to 9.8 with hydrochloric acid.	0.1g 0.2g 97g

Adjust final volume to 1000 ml with distilled water. Store at 4°C

Training Workshop on Diagnostics of Plant Viruses

IPB, University of the Philippines Los Baños | 17-28 August 2015

Total DNA Extraction

STOCK SOLUTIONS <u>1 M Tris-Base</u> Tris-Base Distilled water	60.55 g 500 ml
<u>0.5 M EDTA</u> EDTA Distilled water	18.61 g 100 ml
<u>5 M NaCl</u> Sodium chloride (NaCl) Distilled water	29.22 g 100 ml

Dellaporta Extraction Buffer

Stock Concentration	Final Concentration	100 ml	500 ml
1M Tris-Base pH 8.0	100 mM	10 ml	50 ml
0.5M EDTA	8.5 mM	1.7 ml	8.5 ml
5 M NaCl	500 mM	10 ml	50 ml
Distilled Water	-	78.3 ml	391.5 ml
	Sterilize		
2-mercaptoethanol	10 mM	78 µl	391 µl

5M Potassium Acetate (KAC)	
KAC	49.07 g
Distilled water	100 ml
Sterilize	
20 % SDS	
Sodium dodecyl sulphate (SDS)	10 g
Distilled water	100 ml
Sterilize	

CTAB Extraction Buffer

Stock Concentration	Final Concentration	100 ml	500 ml
1M Tris-HCl pH 8.0	0.1 M	10 ml	50 ml
0.5M EDTA	0.02 M	4 ml	20 ml
5 M NaCl	1.4 M	28 ml	140 ml
СТАВ	2% (w/v)	2 g	10 g
Distilled Water	-	58 ml	290 ml
Sterilize			

<u>TE Buffer</u>

Stock Concentration	Final Concentration	100 ml	1 L
1M Tris-HCl pH 7.5	10 mM	1 ml	10 ml
0.5M EDTA pH 8.0	1 mM	200 µl	2 ml
Distilled Water	-	98.8 ml	988 ml
Sterilize.			

Gel Electrophoresis

<u>10x TI</u> N	BE Buffer Tris Base Boric Acid 0.5 M EDTA Distilled water ote: TAE can also be used.	54.5 g 27.2 g 20 ml 480 ml
<u>0.5x T</u>	<u>BE Buffer</u> 10x TBE buffer Distilled water	50 ml 950 ml
<u>10X L</u>	oading Dye Sucrose 1M Tris-HCl pH 7.5 0.5M EDTA Bromophenol Blue Distilled water	6.5 g 100 μl 200 μl 0.03 g 9.7 ml
<u>6X Lo</u> a	ading Dye 10X Loading dye Distilled water	600 μl 400 μl
<u>1Kb pl</u>	<u>lus DNA Ladder (0.1 μg/ml)</u> DNA Marker (1 μg/ml) 6X Loading Dye Sterile distilled water	50 µl 85 µl 365 µl
<u>1.2% (</u>	<u>Gel</u> Agarose 0.5X TBE buffer Melt agarose using a microwave.	0.48 g 40 ml

Training Workshop on Diagnostics of Plant Viruses

IPB, University of the Philippines Los Baños | 17-28 August 2015

GelRed Solution	
5M NaCI*	2 ml
GelRed	30 µl
Distilled water	98 ml
*optional	

**or add 4 µl GelRed directly to the melted agarose gel (40 ml) Note: Avoid exposing the solution in light as GelRed is light sensitive.

PCR

10 µl
10 µl
10 µl
10 µl
60 µl

Reconstitution of primers

- 1. Centrifuge tubes for a few seconds then add the appropriate volume of TE buffer.
- 2. Rehydrate for 2 min then vortex for 15 sec.
- 3. Short spin to collect contents in the bottom. Store at -20 °C.

Training Workshop on Diagnostics of Plant Viruses

IPB, University of the Philippines Los Baños | 17-28 August 2015

Mechanical Inoculation

Phosphate Buffer Solution A KH ₂ PO ₄ Distilled water Solution B Na ₂ HPO ₄ x 2H ₂ O Distilled water	1.36 g 1000 ml 1.78 g 1000 ml
0.01M Phosphate Buffer, pH 7.2 Solution A Solution B	51 ml 49 ml
Phospahte Buffer (10X) KH ₂ PO ₄ Na ₂ HPO ₄ x 2H ₂ O Distilled water Adjust pH to 7.2 then volume to 1000 ml	2.72 g 14.2 g 800 ml
0.01M Phosphate Buffer, pH 7.2 Phosphate buffer (10X) Distilled water	100 ml 900 ml

ELISA Form



PCR Form for Begomoviruses

	PCR Experimental	Set-up for Bego	omoviruses		
Date:	Ref. Expt. #:		Ref. Gel #		
Title/Description:			Therm	alcycler:	
Time Start:	Time	End:	Total R	un Time:	
	Test	Parameters			
PCR Profile:			DNA Samples:		
Denaturation	94°C1min				
Annealing	57°C 2 min 30x				
Elongation	72°C 2 min				
Final Extension	72°C 10 min ¹				
Primer Pair: PAL1v197	78B/PAR1c715H				
Expected Product Siz	e: ~1.5kb				
Target Gene(;): partie	al DNA-A				
.	PCR COCR				
Components	[)tock]	[Findi]	1X=25µI		
DEPC Water			17.65		
PCR Butter	10x	1x	2.5		
	50mM	2.5mM	1.25		
Drive or E	10mM	0.2mM	0.5		
Primer F	ΙΟμΜ	0.2µM	0.5		
	ΙΟμΜ	0.2µм	0.5		
			0.1		
Total Peaction Volum			2		
	Aliquet 92		25		
Col Elostron horosta	Aliquot 23				
Wei Electrophoresis	Ualtana	Notes/ Kemarks	8		
%del: Dun Timo:	Staining Time				
Ruit fille;	Junning Time:				
	Doc	umentation		· ·	
	200				
Set-up by:					

PCR Form for BBTV

	PCR Experime	ntal Set-up for B	BTV		
Date:	Ref. Expt. #:		Ref. Gel #:		
Title/Description:		Thermalcycler:		ılcycler:	
Time Start:	Time	End:	Total Ru	ın Time:	
	Test	Parameters 199			
PCR Profile:			DNA Samples:		
Initial Denaturation	94°C 4 min	-			
Denaturation	94°C 1 min				
Annealing	61°C 1 min – 29x				
Elongation	72°C 2 min_				
Final Extension	72°C 10 min				
Primer Pair: D11/D12					
Expected Product Siz	e: ~1.1kb				
Target Gene(;): DNA	-R	-			
		<u> </u>			
	·				
	PCR Coch	tail Preparation			
Components	[\$tock]	[Final]	1x=25µI		
DEPC Water			17.65		
PCR Buffer	10x	1x	2.5		
MgCl2	50mM	2.5mM	1.25		
dNTPS	10mM	0.2mM	0.5		
Primer F	10µM	Ο.2 μ Μ	0.5		
Primer R	10µM	0.2µM	0.5		
Таq			0.1		
DNA			2		
Total Reaction Volun	16		25		
	Aliquot 23	to each DNA tube			
Gel Electrophoresis		Notes/Remarks:			
%Gel:	Voltage:				
Run Time:	Staining Time:				
	Doci	umentation			
Set-up by:					

RT-PCR Form for Potyviruses

	PCR Experiment	al Set-up for Po	tyviruses		
Date:	Ref. Expt. #:		Ref. Gel #:		
Title/Description:			Thermalcycler:		
Time Start:	Time	End:	Total R	un Time:	
	Test	Parameters			
PCR Profile:			DNA Samples:		
cDNA Synthesis	55°C 30min				
Initial Denaturation	94°C 2 min				
Denaturation	94°C 15 sec				
Annealing	55°C 30sec – 40x				
Elongation	68°C 1 min _				
Final Extension	68°C 5 min				
Primer Pair: Potyvirid	1/Potyvirid2				
Expected Product Siz	e: ~1.6kb				
Target Gene(\$): 3' term	inal portion of the genome				
	PCR Cocl	ktail Preparation	l		
Components	[\$tock]	[Final]	1x=10µl		
DEPC Water			2.2		
Reaction mix	2x	1x	5		
Primer F	10µM	0.5µM	0.2		
Primer R	10µM	0.5µM	0.2		
Superscript III RT	·	•	0.4		
RNA			2		
Total Reaction Volum	1e		10		
	Aliquot 8	I to each DNA tube	9	1	
Gel Electrophoresis		Notes/Remarks	1		
%Gel	Voltage		•		
Run Time:	Staining Time:				
	Doc	umentation		•	
Set-up by:					

RT-PCR Form for PRSV

PCR Experimental Set-up for PRSV				
Date:	Ref. Expt. #:		Ref. Gel #:	
Title/Description:	-		Thermalcycler:	
Time Start:	Time End:		Total Ru	ın Time:
	Test	Parameters		
PCR Profile:			DNA Samples:	
cDNA Synthesis	55°C 30min			
Initial Denaturation	94°C 2 min			
Denaturation	94°C 15 sec 🗋			
Annealing	58°C 30sec – 40x			
Elongation	68°C 1 min _			
Final Extension	68°C 5 min			
Primer Pair: MB11/MB	12			
Expected Product Siz	e: ~900 bp			
Target Gene(s): CP q	ene			
	PCR Cock	tail Preparation		
Components	[\$tock]	[Final]	1x=10µl	
DEPC Water			2.2	
Reaction mix	2x	1x	5	
Primer F	10µM	0.5µM	0.2	
Primer R	10µM	0.5µM	0.2	
Superscript III RT			0.4	
RNA			2	
Total Reaction Volun	ne		10	
	Aliquot s µl	to each DNA tube		
Gel Electrophoresis		Notes/Remarks:		
%Gel:	Voltage:			
Run Time:	Staining Time:			
	2			
	Doc	umentation	•	
Set-up by:				

DETECTION OF PLANT VIRUSES THROUGH SEROLOGICAL ASSAY (Enzyme-linked Immunosorbent Assay, ELISA)



Enzyme-linked Immunosorbent Assay (ELISA)

 Test that uses antibodies and an enzymemediated color change to identify a substance

- Popular format of 'WET- LAB' type analytic biochemistry assay that uses a solid phase enzyme immunoassay

- Diagnostic tool in medicine, plant pathology as well as quality control check in industries

DETECTION OF PLANT VIRUSES THROUGH SEROLOGICAL ASSAY (Enzyme-linked Immunosorbent Assay, ELISA)





DETECTION OF PLANT VIRUSES THROUGH SEROLOGICAL ASSAY (Enzyme-linked Immunosorbent Assay, ELISA)





DETECTION OF PLANT VIRUSES THROUGH SEROLOGICAL ASSAY (Enzyme-linked Immunosorbent Assay, ELISA)



